

## MATERIAL DATA SHEET

### Tetra-Ub/Ub4 Non-Hydrolyzable Chains (linear), Agarose Cat. # UCN-712

This linear ubiquitin fusion protein is resistant to the activity of enzymes (DUB's) that cleave the peptide linkage between adjacent ubiquitin molecules. Ub is not expressed directly as free Ub but rather as linear fusions either to itself or to certain ribosomal protein subunits. These Ub-fusion precursors are proteolyzed by deubiquitinating enzymes (DUBs) at the appropriate junction points to yield active Ub monomers with C-termini ending in GG. There are likely several intracellular DUBs which perform this essential processing role. This product may be useful in analyzing interactions between linear ubiquitin and proteins that contain ubiquitin-associated domains (UBAs) or ubiquitin-interacting motifs (UIMs).

#### Product Information

<b>Quantity:</b>	250 µl resin supplied as a 50% slurry in 20% ethanol.
<b>Stock:</b>	Non-hydrolyzable Ub4 (linear) chains coupled to beads at 2 mg/ml (58 nmol/ml).

#### Use and Storage

<b>Use:</b>	Prepare resin by rinsing with 10 volumes of water to remove ethanol storage buffer. Equilibrate resin by washing with 10 volumes of desired start buffer. Binding and elution of material is dependent on individual experimental conditions and requirements.
<b>Storage:</b>	Polyubiquitin-agarose can be re-used if properly maintained. After use, clean resin with a wash cycle of 5 volumes 100 mM HEPES pH 8.0, 500 mM NaCl followed by 5 volumes 100 mM NaOAc pH 4.5, 500 mM NaCl. Repeat twice, then rinse resin with a low salt buffer. Store resin at 4°C in neutral aqueous buffer containing 1 mM NaN <sub>3</sub> or 20 % ethanol as a preservative. DO NOT FREEZE.

#### Literature

<b>References:</b>	Gerlach B., <i>et al</i> (2011) <u>Nature</u> <b>471</b> : 591-596 Emmerich C., <i>et al</i> (2011) <u>Science Signaling</u> <b>4</b> : re5 Sato Y., <i>et al</i> (2011) <u>PNAS</u> <b>108</b> : 20520-20525
--------------------	---

*For Laboratory Research Use Only, Not For Use in Human*

840 Memorial Drive, Cambridge, MA 02139 Phone: 617-576-2210 FAX: 617-492-3565

[www.bostonbiochem.com](http://www.bostonbiochem.com)

The contents of this datasheet (unless otherwise noted) are Copyright © 2012 Boston Biochem, Inc. All rights reserved. Duplication in whole or in part is strictly prohibited without the express written consent of Boston Biochem, Inc. "Boston Biochem" is a Trademark of Boston Biochem, Inc., 840 Memorial Drive, Cambridge, MA 02139.