Lot # XXXXX



An R&D Systems Company

MATERIAL DATA SHEET

Poly-SUMO2 Hydrolysis-Resistant Chains (3-8) Cat. # ULN-220

Sentrin protease (SENP)-resistant poly-SUMO2 chains can be used to investigate mechanisms of chain recognition by SUMO-targeted ubiquitin ligases (STUbl's) or other proteins that contain SUMO-interacting domains (SIM's). These K11-linked chains are formed enzymatically with recombinant SUMO2 containing a mutation of glutamine 90 to proline, which renders them approximately 500-fold more resistant to disassembly by SENP's than are wild-type chains. Mono- and di-SUMO2 have been removed from the chain mixture.

Product Information

Quantity: 50 μg

Stock: X mg/ml in 50 mM Hepes pH 8.0, 100 mM NaCl, 1 mM DTT.

Purity: > 90% by SDS-PAGE

Use & Storage

Use: Typical concentrations will depend on specific assay conditions and method of

detection.

Storage: Store at -20° C. Avoid multiple freeze/thaw cycles.

Literature

References: Bekes, M., et al. (2011), J. Biol. Chem. **286**: 10238-10247

Bylebyl G.R., et al. (2003) <u>J. Biol. Chem.</u> **278**:44113-44120 Dohmen R. J. (2004) <u>Biochem. Biophys. Acta.</u> **1695**:113-131

Johnson E. (2004) Ann. Rev. Biochem. 73:355-382

Tatham M.H., et al. (2001) J. Biol. Chem. 276: 35368-35374

For Laboratory Research Use Only, Not For Use in Humans

Rev: 07/19/2012

840 Memorial Drive, Cambridge, MA 02139 Phone: 617-241-7072 FAX: 617-492-3565 www.bostonbiochem.com