

Lot # XXXXX

BostonBiochem

MATERIAL DATA SHEET

Human SUMO3, Monoclonal Antibody Cat. # A-720

Small Ubiquitin-like Modifiers (SUMOs) are a family of small, related proteins that can be enzymatically attached to a target protein by a post-translational modification process termed SUMOylation. Unlike ubiquitination, which targets proteins for degradation, SUMOylation participates in a number of cellular processes, such as nuclear transport, transcriptional regulation, apoptosis, and protein stability. All human SUMO proteins share a conserved ubiquitin-like domain and a C-terminal diglycine cleavage/attachment site. Human SUMO3, also known as SMT3A, is synthesized as a 103 amino acid (aa), 11.5 kDa propeptide that contains an 11 aa C-terminal prosegment. Following prosegment cleavage, the C-terminal glycine may be enzymatically attached to a lysine on a target protein. Human SUMO3 shares 83% sequence identity with SUMO3 from mouse. SUMO3 also has very high sequence homology to SUMO2 and SUMO4, 87% and 75%, respectively. SUMO3 shares only 47% sequence identity to SUMO1.

Product Information

Quantity:	100 µg
Source:	Monoclonal Rat IgG _{2A} Clone # 907614
Antigen:	Peptide sequence from SUMO3. Accession Number P55854
Purification:	Protein G purified from hybridoma culture supernatant
Stock:	0.5 mg/mL in PBS, pH 7.4, 50% glycerol, 0.09% sodium azide

Use & Storage


Specificity:	This antibody detects endogenous, human SUMOylated proteins in Western blots. This antibody detects purified, recombinant SUMO3, but has no cross-reactivity to recombinant human SUMO1, 2, or 4.
Use:	Recommended concentration for Western Blot is 0.1- 0.5 µg/ml.
Storage:	Store at -20°C.

840 Memorial Drive, Cambridge, MA 02139 Phone: 617-576-2210 FAX: 617-492-3565

www.bostonbiochem.com

The contents of this datasheet (unless otherwise noted) are Copyright © 2008 Boston Biochem, Inc. All rights reserved. Duplication in whole or in part is strictly prohibited without the express written consent of Boston Biochem, Inc. "Boston Biochem" is a Trademark of Boston Biochem, Inc., 840 Memorial Drive, Cambridge, MA 02139.

Sample Western Blot Data



10-40 ng of recombinant pro-SUMO proteins were separated on a reducing 10-20% SDS-PAGE gel. Western blots were developed using PVDF membranes and α -SUMO3 (A-720) mAb primary at 0.5 μ g/ml followed by HRP-labeled anti-rat (R&D Systems # HAF008) secondary antibody at 1:2000 dilution. Only the SUMO3 protein was detected in the blot.

Literature

- References:** Desterro J.M., *et al.* (1997) FEBS. Lett. **417**: 297-300
 Okama T., *et al.* (1999) Biochem. Biophys. Res. Comm. **254**: 693-698
 Seeler J-S. and Dejean A. (2003) Nat. Rev. **4**: 690-699
 Su H-L., *et al.* (2002) Gene **296**: 65-73
 Tatham M.H., *et al.* (2001) J. Biol. Chem. **276**: 35368-35374
 Yeh E.T.H., *et al.* (2000) Gene **248**: 1-14

For Laboratory Research Use Only, Not For Use in Humans

Rev: 02/14/2014

840 Memorial Drive, Cambridge, MA 02139 Phone: 617-576-2210 FAX: 617-492-3565

www.bostonbiochem.com

The contents of this datasheet (unless otherwise noted) are Copyright © 2008 Boston Biochem, Inc. All rights reserved. Duplication in whole or in part is strictly prohibited without the express written consent of Boston Biochem, Inc. "Boston Biochem" is a Trademark of Boston Biochem, Inc., 840 Memorial Drive, Cambridge, MA 02139.