

Lot # XXXXX

BostonBiochem

## MATERIAL DATA SHEET

### 10X Ubiquitin Reconstitution Buffer

Cat. # B-90

This buffer is formulated to provide the optimal formulation for the reconstitution and subsequent dilution of lyophilized Ubiquitin and Ubiquitin-like proteins.

#### Product Information

**Quantity:** 5 X 1 ml vials

**Stock:** 10X Solution

#### Use & Storage

**Use:** Dilute to 1X in deionized, sterile water prior to addition to lyophilized ubiquitin. We strongly recommend that lyophilized materials be centrifuged at 14k xg for 2-5 minutes prior to reconstituting to pellet all material to the bottom of the tube.

**Storage:** Store at -20°C. Mix well prior to use.

*For Laboratory Research Use Only, Not For Use in Humans*

Rev: 01/22/2015

840 Memorial Drive, Cambridge, MA 02139 Phone: 617-576-2210 FAX: 617-492-3565

[www.bostonbiochem.com](http://www.bostonbiochem.com)

The contents of this datasheet (unless otherwise noted) are Copyright © 2008 Boston Biochem, Inc. All rights reserved. Duplication in whole or in part is strictly prohibited without the express written consent of Boston Biochem, Inc. "Boston Biochem" is a Trademark of Boston Biochem, Inc., 840 Memorial Drive, Cambridge, MA 02139.