

**MATERIAL DATA SHEET****Ramos Cell 20S Immunoproteasome Protein, *human***  
**Cat. # E-371**

The 20S Immunoproteasome is a modified form of the constitutively active 20S Proteasome core particle and is the catalytic subunit of the multi-complex Immunoproteasome. The structure of the 20S Immunoproteasome is similar to the 20S Proteasome, which is composed of 28 non-identical subunits arranged into four stacked rings. However, during 20S Immunoproteasome assembly, the three catalytic beta subunits,  $\beta 1$ ,  $\beta 2$ , and  $\beta 5$ , in the two interior rings of the 20S Proteasome are replaced by three IFN $\gamma$ -inducible catalytic subunits:  $\beta 1i$ /LMP2,  $\beta 2i$ /MECI-1, and  $\beta 5i$ /LMP7. The 20S Immunoproteasome is commonly associated with the 19S, PA28  $\alpha/\beta$ , or the PA28 $\gamma$  regulatory complexes. 20S Immunoproteasome expression is enriched in antigen presenting cells of the immune system where the 20S Immunoproteasome selectively degrades intracellular proteins in a manner that optimizes the generation of peptides for MHC class I antigen presentation. Selective inhibition of 20S Immunoproteasome proteolytic activity using small molecule inhibitors is being examined for therapeutic intervention in cancer and inflammatory diseases. The protein complex has been purified from Burkitt's lymphoma cells (Ramos, RA1, CRL-1596) following stimulation with IFN- $\gamma$  and IL2.

**Product Information**

<b>Quantity:</b>	25 $\mu$ g
<b>Stock:</b>	X mg/ml (X $\mu$ M) in 50 mM HEPES pH 7.6, 100 mM NaCl, 1 mM DTT.
<b>MW:</b>	700 kDa
<b>Purity:</b>	> 95% by SDS-PAGE under reducing conditions and visualized by Colloidal Coomassie Blue stain.

**Use & Storage**

<b>Use:</b>	The human 20S Immunoproteasome is able to degrade substrates in an ATP-independent manner. It can be activated chemically with SDS (0.035%) or by the addition of PA28. Reaction conditions will need to be optimized for each specific application. We recommend an initial 20S Proteasome concentration of 0.5-5 nM.
<b>Storage:</b>	Store at -80 $^{\circ}$ C. Avoid multiple freeze/thaw cycles.

## Literature

- References:** Cascio, P. *et al.* (2001) EMBO J. **20**: 2357  
Ferrington, D.A. & D.S. Gregerson (2012) Prog. Mol. Biol. Transl. Sci. **109**: 75  
Kim, H.M. *et al.* (2011) Biochim. Biophys. Acta **1809**: 67  
Kloetzel, P.M. (2001) Nat. Rev. Mol. Cell Biol. **2**: 179  
Lee, W. & K.B. Kim (2011) Curr. Top. Med. Chem. **11**: 2923  
Niewerth, D., *et al.* (2014) J Hema & Onc. **7**:7  
Stadtmueller, B.M. & C.P. Hill (2011) Mol. Cell **41**: 8  
Xie, Y. (2010) J. Mol. Cell Biol. **2**: 308□

*For Laboratory Research Use Only, Not For Use in Humans*

Rev: 04/8/2016

840 Memorial Drive, Cambridge, MA 02139 Phone: 617-576-2210 FAX: 617-492-3565  
[www.bostonbiochem.com](http://www.bostonbiochem.com)

The contents of this datasheet (unless otherwise noted) are Copyright © 2008 Boston Biochem, Inc. All rights reserved. Duplication in whole or in part is strictly prohibited without the express written consent of Boston Biochem, Inc. "Boston Biochem" is a Trademark of Boston Biochem, Inc., 840 Memorial Drive, Cambridge, MA 02139.