

**MATERIAL DATA SHEET****FAT10 Agarose, human recombinant**  
**Cat. # UL-920**

FAT10 covalently coupled to agarose beads via primary amines allowing for a fully functional C-terminus. Useful for isolation and capture of FAT10 interacting proteins and/or enzymes that have an affinity for this ubiquitin-like protein.

**Product Information**

<b>Quantity:</b>	0.5 ml
<b>Stock:</b>	0.5 ml resin supplied in a 1 ml total volume of 50 mM Hepes pH 7.5, 250 mM NaCl.

**Use & Storage**

<b>Use:</b>	Equilibrate resin by washing with 5-10 ml desired start buffer. Binding and elution of material is dependent on individual experimental conditions.
<b>Storage:</b>	The agarose can be re-used for at least 5-10 applications if properly maintained. After use, clean resin with 5ml 50 mM Tris pH 9.0, 1M KCL. Remove cleaning solution by washing resin with 5 ml storage buffer. Resin should be stored at 4°C and 0.01% sodium azide can be used as a bacteriostatic agent. DO NOT FREEZE.

**Literature**

<b>References:</b>	Bates E.E. <i>et al.</i> (1997) <u>J. Immunol.</u> <b>27</b> : 2471-2477
	Canaan A. <i>et al.</i> (2006) <u>Mol.Cell.Biol</u> <b>27</b> : 5180-5189
	Fan W. <i>et al.</i> (1996) <u>Immunogen.</u> <b>44</b> : 97-103
	Hipp M.S. <i>et al.</i> (2004) <u>J. Biol Chem.</u> <b>279</b> : 16503-16510
	Hipp M.S., <i>et al.</i> (2005) <u>Mol. Cell. Biol.</u> <b>25</b> : 3483-3491
	Lee C.G. <i>et al.</i> (2003) <u>Oncogene.</u> <b>22</b> : 2592-2603
	Liu Y.C. <i>et al.</i> (1999) <u>Proc. Natl. Acad. Sci.</u> <b>96</b> : 4313-4318
	Lim C.B. <i>et al.</i> (2006) <u>Cell. Div.</u> <b>1</b> :20
	Raasi S. <i>et al.</i> (2001) <u>J. Biol Chem.</u> <b>276</b> : 35334-35343
	Ross M.J. <i>et al.</i> (2006) <u>.Am.Soc.Nephrol.</u> <b>17</b> : 996-1004

***For Laboratory Research Use Only, Not For Use in Humans***