
MATERIAL DATA SHEET

**Poly-SUMO3 Wild-type Chains (2-8)
Cat. # ULC-310**

Poly-SUMO-3 chains can be used to investigate mechanisms of chain recognition, binding and hydrolysis by SUMO-specific isopeptidases (SENPs), SUMO-specific E3 ligases or other proteins that contain SUMO-3 binding domains. This product is formed enzymatically with wild-type human recombinant SUMO-3 linked via lysine 11 which is the point of attachment for the C-terminal glycine of the preceding SUMO-3. Mono-SUMO-3 has been removed from the chain mixture.

Product Information

Quantity:	25 µg
Stock:	X mg/ml in 50 mM Hepes pH 8.0, 100 mM NaCl, 1 mM DTT. Concentration will vary with specific Lot #.
Purity:	> 95% by SDS-PAGE

Use & Storage

Use:	Typical concentrations will depend on specific assay conditions and method of detection.
Storage:	Store at -80°C. Avoid multiple freeze/thaw cycles.

Literature

References:	Bylebyl G.R., <i>et al.</i> (2003) <i>J. Biol. Chem.</i> 278 :44113-44120 Dohmen R. J. (2004) <i>Biochem. Biophys. Acta.</i> 1695 :113-131 Johnson E. (2004) <i>Ann. Rev. Biochem.</i> 73 :355-382 Tatham M.H., <i>et al.</i> (2001) <i>J. Biol. Chem.</i> 276 : 35368-35374
--------------------	---

For Laboratory Research Use Only, Not For Use in Humans